




OtoVet

HIGH PERFORMANCE VIDEO OTOSCOPE

Take the best video otoscopy system into your next exam. BioVision Veterinary Endoscopy's OtoVet Video Otoscopy system gives you the outstanding image performance of the **SurgVIEW** Xenon Imaging System, an otoscope with a 90mm working length and a great picture.



Elbows | Shoulders | Knees | Hips | Ultra Portable | Fully Integrated MediaCapture*



OtoVet

HIGH PERFORMANCE VIDEO OTOSCOPE

OtoVet High Performance Video Otoscope Part No. BVTSUITE-00005

| Item | Description |
|-----------|--|
| SurgView | Xenon imaging system w/ MediaCapture USB |
| Camera | High resolution (480 Line) w/ integrated fiber optic cable |
| Coupler | 20mm F.L. C-mount coupler |
| Otoscope | 0 deg. 5.5-20mm O.D., 90mm W.L., 2.1mm channel |
| Forceps | 5F. flexible alligator grasping 280mm W.L. |
| Forceps | 5F. flexible biopsy, 280mm W.L. |
| Curette | Flexible, 213mm W.L. |
| Accessory | Cleaning brush for working channel |
| Training | Training |

SurgVIEW

SurgASSIST

MediaCaptureUSB

BIO
VISION Technologies, LLC
VETERINARY
ENDOSCOPY

BVT93062-00001-G
062409