

Less

(bulk and complication)


is

More.

(portable and user-friendly)

Introducing BioVision Veterinary Endoscopy's EndoDiagnostic and Surgical Suite.

From general video otoscopy to advanced laparoscopy, the EDSS scales nicely to cover many clinical needs for rigid endoscopy. Not only do you get proven tools, but you also get proven techniques to make rigid endoscopy a part of daily practice.



EDSS

EndoDiagnostic+Surgical

S U I T E



EDSS

EndoDiagnostic+Surgical
S U I T E

EndoDiagnostic and Surgical Suite Part No. BVTSUITE-00001

Item

SurgView
Camera
SurgAssist
Laparoscope
Otoscope
Coupler
EndoBlade (2x)
Trocar (2x)
Trocar
Forceps
Forceps
Handpiece
Accessory
Hose
Footswitch
Regulator
Hose
Plate
Wires
Accessory
Accessory
DVD
Training

Description

Xenon imaging system w/ MediaCapture USB
High resolution (480 Line) w/ integrated fiber optic cable
Insufflator/electrosurgical generator
11mm O.D. 6mm I.D., 270mm W.L. 45 deg oculus
0 deg. 5.5-20mm O.D., 90mm W.L., 2.1mm channel
20mm F.L. C-mount coupler
5.16mm O.D. , 420 W.L., Bi-polar cautery instrument
11/12mm, 100mm W.L.
5mm, 75mm W.L.
5mm O.D., 430mm W.L. babcock w/ ratchet
5mm O.D., 430mm W.L. cup biopsy w/ tooth
Monopolar pencil for extracorporeal use
Monopolar connecting wire for endoscopic instruments
Patient insufflation, open end to male luer lock
Electrosurgical generator activator switch
E-tank saddle w/ pressure gauge
CO2 supply, regulator to BV quick-connect
Grounding plate for monopolar cautery
Wires for grounding plate
Cleaning brush for operating laparoscope
Cleaning brush for otoscope
Procedure DVD & client educational materials
Hands- on surgical training

Options

Table "TippyTable" surgical positioning table
Forceps 5mm O.D., 330mm W.L. babcock w/ ratchet
Forceps 10mm O.D., 330mm W.L. babcock w/ ratchet
Sterile Tray Clear sterilization tray, 23" x 8" x 3.75"
Sterile Drapes Sterile drapes for camera w/ operating laparoscope
Cart/Monitor/Arm Adjustable cart w/ AC power & swivel arm w/ 19" LCD display

Surgical Training

Hands-on surgical endoscopic training is included with each EndoDiagnostic & Surgical Suite. Some of the most progressive surgical techniques are taught in 8 hours of intensive training.

SurgVIEW™

SurgASSIST™

MediaCapture^{USB}

BIOVISION Technologies, LLC
VETERINARY
ENDOSCOPY

BVT93062-00003-G
062409