

BioVision Veterinary Endoscopy

Rigid Operating Endoscopes

Custom built, high quality rigid endoscopes for diagnostic and therapeutic procedures. You will find these highly reliable scopes becoming an essential part of your daily practice.



	O.D.	I.D.	W.L.	Oblique	Max. Inst.	Part Number
Veterinary Otoscope	5.4-19.97mm	1.9mm	90mm	0°	5 French	BVTENDO3-00001
MultiPurpose Rigid	2.7mm	-	270mm	30°	-	BVTENDO3-00007
Diagnostic Sheath	5.5mm	-	224mm	-	-	BVTINST2-00005
Operating Sheath	6.1mm	2mm	190mm	-	5 French	BVTINST2-00006
Obturator	-	-	-	-	-	BVTINST3-00005
Operating Laparoscope	11mm	6mm	270mm	0°	5.8mm	BVTENDO1-00001



BioVision Veterinary Endoscopy, 221 Corporate Circle Suite H, Golden, Colorado 80401

biovisiontech.com

BVT93062-00008-B
040605