




# Arthro<sup>+</sup> *mobility Advance* **Veterinary Arthroscopy** E N D O S C O P Y   S Y S T E M

Think of it: Your whole arthroscopy system in one easy to roll case. It's all in there, SurgView Xenon Imaging System, sterile tray, cords, and a cubby for sterile hardware packs. Pop the handle up and roll it into your next appointment. The reinforced, hardcase is nearly indestructible and custom cut foam holds all the contents securely. Yes, now you can take it with you.



**Arthro<sup>+</sup>**  
*mobility Advance*  
**Veterinary Arthroscopy**  
ENDOSCOPY SYSTEM

**ArthroPlus Mobility Advance Veterinary Arthroscopy Endoscopy System**  
**Part No. BVTSUITE-00010**

Item	Part Number	Description
Case	BVTACCE1-00003	Black resin fiber with casters and telescoping handle
Accessory	BVTACCE1-00004	Custom Foam for SurgView, Sterile Tray and Cords/etc.
Accessory	BVTACCE1-00011	Custom Foam for Arthroscopy Case
Tray	BVTACCE1-00012	Bi-level sterilization with silicone mat and instrument holders, 12.5" x 10" x 5.75
SurgView	BVT01162-01001	Xenon Imaging System with MediaCapture USB
Camera	BVT01163-01001	High Resolution (480 Line)/Fiber Optic Cable Combination
Coupler	BVTACCE1-00001	28mm Imaging EndoCoupler with locking spin and quick release
Arthroscope	BVTENDO3-00006	30 deg, 1.9mm O.D., 60mm W.L.
Cannula	BVTINST2-00004	3.97mm O.D., 50mm W.L.
Obturator	BVTINST3-00006	Conical Blunt**
Arthroscope	BVTENDO3-00003	30 deg, 2.7mm O.D., 140mm W.L.
Cannula	BVTINST2-00002	2.8mm O.D., 45mm W.L.
Obturator	BVTINST3-00001	Conical Blunt**
Forceps	BVTINST5-00011	3.5mm, Meniscus grasping with ratchet
Forceps	BVTINST5-00010	2.0mm, Rigid grasping forcep
Curette	BVTINST5-00007	2.0mm Spoon
Probe	BVTINST5-00009	90 deg, hook style
Handle	BVTINST5-00008	Handle for hook probe
Accessory	BVTACCE1-00016	Cleaning brush for working channel
Training	BVTACCE1-TRAIN	On-Site Training

\*\* Other styles available: Pyramidal Sharp, Spherical Blunt

