


# Arthro+ Advance Veterinary Arthroscopy

E N D O S C O P Y   S Y S T E M

## Quality, Reliability, Performance.

These are the cornerstones of your practice and the cornerstones on which we built BioVision Veterinary Endoscopy's ArthroPlus Veterinary Arthroscopy System. The two ultra high quality arthroscopes take advantage of our SurgView Xenon Imaging System's superior light and camera imaging systems. A set of essential hand tools allow you to get to work without having to piece a system together.

Our MediaCapture USB\* system is a fully-integrated media solution that allows you to capture stills, audio tagged stills and video from the sterile field. Additionally, you can stream video over a USB connection for live capture.



Elbows | Shoulders | Knees | Hips | Ultra Portable | Fully Integrated MediaCapture\*

# Arthro<sup>+</sup> Advance Veterinary Arthroscopy

ENDOSCOPY SYSTEM



## ArthroPlus Advance Veterinary Arthroscopy Endoscopy System

Part No. BVTSUITE-00007

Item	Part Number	Description
SurgView	BVT01162-01001	Xenon Imaging System with MediaCapture USB
Camera	BVT01163-01001	High Resolution (480 Line)/Fiber Optic Cable Combination
Coupler	BVTACCE1-00001	28mm Imaging EndoCoupler with locking spin and quick release
Arthroscope	BVTEND03-00006	30 deg, 1.9mm O.D., 60mm W.L.
Cannula	BVTINST2-00004	3.97mm O.D., 50mm W.L.
Obturator	BVTINST3-00006	Conical Blunt**
Arthroscope	BVTEND03-00003	30 deg, 2.7mm O.D., 140mm W.L.
Cannula	BVTINST2-00002	2.8mm O.D., 45mm W.L.
Obturator	BVTINST3-00001	Conical Blunt**
Forceps	BVTINST5-00011	3.5mm, Meniscus grasping with ratchet
Forceps	BVTINST5-00010	2.0mm, Rigid grasping forceps
Curette	BVTINST5-00007	2.0mm Spoon
Probe	BVTINST5-00009	90 deg, hook style
Handle	BVTINST5-00008	Handle for hook probe
Accessory	BVTACCE1-00016	Cleaning brush for working channel
Training	BVTACCE1-TRAIN	On-Site Training

\*\* Other styles available: Pyramidal Sharp, Spherical Blunt

**SurgVIEW**™ MediaCapture<sup>USB</sup>