

# Less

(bulk and complication)

# is

# More.

(portable and user-friendly)



**EDSS**  
EndoDiagnostic+Surgical  
S U I T E

## Introducing BioVision Veterinary Endoscopy's EndoDiagnostic and Surgical Suite.

From general video otoscopy to advanced laparoscopy, the EDSS scales nicely to cover many clinical needs for rigid endoscopy. Not only do you get proven tools, but you also get proven techniques to make rigid endoscopy a part of daily practice.



# EDSS

EndoDiagnostic+Surgical  
S U I T E

## EndoDiagnostic and Surgical Suite

Part No. BVTsuite-00001

### Suite Contents:

Item	Description
SurgView	Xenon Imaging System with MediaCapture USB (Available Q1 2005)
Camera Head for SurgView	High Res. Camera with hybrid cable
SurgAssist	Insufflator/Electrosurgical Generator
Operating Laparoscope	11mm O.D. 6mm I.D., 270mm W.L. 45 deg oculus
Otoscope	0 deg. 5.4-25mm O.D., 90mm W.L., 2.1mm channel
Coupler	28mm F.L. C-Mount to B-Style
EndoBlade	5.16mm O.D. , 420 W.L.
Trocar	10/11mm, 100mm W.L.
Trocar	5mm, 75mm W.L.
Trocar Sleeve	18mm, with switching stick
Forceps	5mm endoscopic babcock
Forceps	10mm endoscopic babcock
Forceps	5mm endoscopic cup biopsy with tooth
Footswitch	Electrosurgical generator activator switch
Regulator	E-tank saddle with pressure gauge
Hose	CO <sub>2</sub> Supply, regulator to BV QuickConnect
Hose	Patient insufflation, open end to male luer lock
Training	8 hours on-site training

## Surgical Training

We come to you. On-site surgical endoscopic training is included with each EndoDiagnostic and Surgical Suite. Some of the most progressive surgical techniques are taught in 8 hours of intensive training geared toward both DVMs and their staff.

**SurgVIEW™** **SurgASSIST™** **MediaCaptureUSB**