



EDSS *mobility*

EndoDiagnostic+Surgical

S U I T E

Now the unprecedented compact size of BioVision Veterinary Endoscopy's SurgView and SurgAssist and the whole EndoDiagnostic and Surgical Suite is on the go. The EDSS Mobility is the answer for doctors needing a surgical suite in the field or experiencing increasing demand in their mobile practice. Two virtually indestructible cases hold a whole EDS Suite including all the cords and extras. Surgical endoscopy is now truly portable.

EndoDiagnostic and Surgical Suite - Mobility
Part No. BVTSUITE-00002

Case 1

Item	Part Number	Description
Case	BVTACCE1-00003	Black Resin Fiber with casters and telescoping handle
Accessory	BVTACCE1-00004	Custom Foam for SurgView, Sterile Tray and Cords/etc.
Tray	BVTACCE1-00006	Bi-level sterilization with Si mat and inst. holders, 28.5" x 11.25" x 7.5"
SurgView	BVT01162-01001	Xenon Imaging System with MediaCapture USB
Camera	BVT01163-01001	High Res. Camera with hybrid cable
Laparoscope	BVTENDO1-00001	11mm O.D. 6mm I.D., 270mm W.L. 45 deg oculus
Otoscope	BVTENDO3-00001	0 deg. 5.4-25mm O.D., 90mm W.L., 1.9mm channel
Coupler	BVTACCE1-00026	30 mm F.L. C-Mount to B-Style
EndoBlade (2x)	BVTINST4-00001	5.16mm O.D., 420 W.L.
Trocar (2x)	BVTINST1-00001	10/11mm, 100mm W.L.
Trocar	BVTINST1-00002	5mm, 75mm W.L.
Forceps	BVTINST5-00001	5mm endoscopic babcock
Forceps	BVTINST5-00002	10mm endoscopic babcock
Forceps	BVTINST5-00003	5mm endoscopic cup biopsy with tooth
Handpiece	BVTINST4-00002	Monopolar Pencil for extracorporeal use.
Accessory	BVTINST4-00003	Monopolar connecting wire for endoscopic instruments
Hose	BVTACCE1-00009	Patient insufflation, open end to male luer lock

Case 2

Case	BVTACCE1-00003	Black Resin Fiber with casters and telescoping handle
Accessory	BVTACCE1-00005	Custom Foam for SurgAssist, Regulator, Footswitch, Hoses, Cords/etc.
SurgAssist	BVT01165-01001	Insufflator/Electrosurgical Generator
Footswitch	BVT22800-05111	Electrosurgical generator activator switch
Regulator	BVT65400-01940	E-tank saddle with pressure gauge
Hose	BVT65765-27201	CO2 Supply, regulator to BV QuickConnect
Plate	BVT45065-45050	Grounding plate for monopolar cautery
Wires	BVTACCE1-00008	Wires for grounding plate
Accessory	BVTACCE1-00025	Cleaning brush for operating laparoscope
DVD	BVTVIDEO-00001	Marketing Starter Kit
Training	BVTACCE1-TRAIN	8 hours on-site surgical training

Surgical Training

We come to you. On-site surgical endoscopic training is included with each EndoDiagnostic and Surgical Suite. Some of the most progressive surgical techniques are taught in 8 hours of intensive training geared toward both DVMs and their staff.

SurgVIEW™ **SurgASSIST™** **MediaCaptureUSB**