


Rigid Endoscopes

Custom built, high quality rigid endoscopes for diagnostic and therapeutic procedures.

You will find these highly reliable scopes becoming an essential part of your daily practice.



	O.D.	W.L.	Oblique	Feature
Arthroscope / Feline Rhinoscope	1.9mm	60.5mm	30°	
Cannula	2.78mm	36.7mm	-	1 Stopcock N/R
Arthroscope / Rhinoscope	2.7mm	75mm	30°	
Cannula	3.97mm	50mm	-	1 Stopcock N/R
Arthroscope / Rhinoscope	2.7mm	140mm	30°	
Cannula	3.97mm	105		2 Stopcock Rot.
MultiPurpose Rigid	2.7mm	270mm	30°	
Diagnostic Sheath	5.5mm	224mm	-	2 Stopcock Rot.
Operating Sheath	6.1mm	190mm	-	2 SC, Rot., 2mm chan.
Arthroscope	3mm	150mm	30°	
Cannula	4.4mm	114mm	-	2 Stopcock Rot.
Arthroscope	4mm	185mm	30°	
Cannula	5.8mm	140mm	-	2 Stopcock Rot.
Laparoscopes	5mm	290mm	0°	
	5mm	290mm	30°	
	10mm	320mm	0°	
	10mm	320mm	30°	

Obturator available for all cannulae and sheaths.

